



## PRODUCT DATA SHEET AND USER GUIDE

### Vascupaint™ Silicone Rubber Injection Compounds

#### Description

Vascupaint™ lead-free silicone injection compounds fill and opacify the microvasculature structures for 3D ex vivo imaging using micro-CT, dissection or brightfield microscopes, or SWIR fluorescence imaging.

The injection compounds are available in two different colors including yellow (SKU# MDL-121) for micro-CT imaging and green (SKU# MDL-122) for NIR or SWIR fluorescence imaging as well as clear (SKU# MDL-123) for your custom preparations. Other colors are available upon request.

Vascupaint™ is a less viscous alternative to conventionally used lead-based injection compounds making it easier to perfuse non-surviving organisms as well as human organs for ex vivo 3D vascular cast analysis.

#### Features & Benefits

- X-ray attenuation of bismuth nanoparticles in Vascupaint yellow product allows for high resolution micro-CT imaging
- The sub 750 nm particle size of the bismuth particles in yellow product allows for adequate filling of capillaries and microvasculature
- Low viscosity of mixture allows for perfusion into non-surviving organism at physiological pressures and complete filling of microvasculature
- The catalyst is a defined media that uses water molecules to drive cross-linking of monomers. The bottle of catalyst does not require shaking prior to use and there is no settling of cross-linker and catalyst and the catalyst bottle does not require shaking prior to use.

#### Applications

- Vascular corrosion casting, imaging of angiogenesis, quantification of vascular volumes in ischemia and tumors, microCT analysis of

vascular morphometry, assessing vascular architecture of organs, teaching adjunct, visualization of microcirculation for surgery techniques

#### Vascupaint Kit Components

- 1 bottle of silicone (Green, Yellow, Lime, Blue, Red or Clear), 200 ml
- 1 bottle of diluent, 200 ml
- 1 bottle of catalyst, 20 ml

#### Storage

All components of this kit should be stored at room temperature and are stable for one year from date of manufacture

#### Product Safety and Handling

This product is for R&D use only, not for drug, household, or other uses. Please review the material safety datasheet (MSDS) for proper safety and handling procedures.

#### Preparing Stock Solutions of Silicone and Diluent Mixtures

1. Add 20 ml of Vascupaint to a 50 ml conical tube. Do not use a pipette for this step.
2. Add 25 ml diluent to the same tube
3. Gently invert and store for future use

#### Using Previously Prepared Stock Solutions of Silicone and Diluent Mixtures

Some 'soft settling' of pigments will occur in the stock solution. You can remove this by gently inverting the bottle several times. The pigments will remain suspended in solutions for several hours after this gentle inversion step.

#### Protocol Example

A recommended starting protocol for mouse perfusions is as follows:

1. Gently invert and/or shake the stock solution bottle.
2. Add 9 ml of the stock solution to a 14 ml conical tube or to petri dish for mixing 0.45 ml of catalyst.
3. Cut the needle tip to make it blunt before insertion into the heart. Insert the syringe needle into heart through the left ventricle and all the way up to the aorta. Cut open the right atrium to allow blood and contrast agents to flow out of circulation. It is also possible to cut the edges of the liver instead of the atrium to ensure that the outflow of the contrast occurs less rapidly. Some users have found that clamping away the non-needed branches (e.g. abdominal aorta below renal arteries in case of kidney perfusion or mesenteric vessels for perfusion of the hind limbs or the left ventricle for the brain) increases the quality of the perfusion. See example video of this procedure. Ref. <https://youtu.be/lcm1OG-vUI4>
4. If IP injection of heparin was not performed in step 2, Perfuse with heparinized PBS (5ml) first, and then switch to a syringe, which has been filled with formalin (5 ml) for fixation.
5. Add 0.45 ml (5%) catalyst to the petri dish with stock solution and aspirate mix several times with pipette. Transfer the mixture to a syringe.
6. Following fixation, switch to the syringe with the vascupaint mixture. Before injecting Vascupaint, put a hemostat to clip on the heart to hold the needle in place and add crazy glue to make sure the clip is held in place and that a back-pressure is sensed when introducing the mixture into the organism. By doing so, the needle will be fixed in the aorta and leakage of Vascupaint mixture out of the heart will be minimized, if not completely avoided. A perfusion pump can be used to administer Vascupaint mixture slowly to avoid any leakages or disruptions of the hemostat clamp. Note: avoid creating air bubbles at the tip of the needle while switching syringes.
7. Leave the mice in the fridge overnight. For brain perfusion procedures, leave the mice upside down at 4°C, overnight.
8. The following day, the Vascupaint inside the mouse's body should be hardened and organs can be removed carefully. For brain procedures, the head and exposing the skull. The organ can be further processed using standard tissue digestion techniques for analysis with dissecting or brightfield microscopy (\*this application may require an increased level of catalyst to ensure proper hardening of the Vascupaint). The intact organ is ready for analysis with microCT imaging for the Yellow version of Vascupaint and surface vessels can be assessed with visual inspection or microscopy to determine success of perfusion.