

SOMNI  
SCIENTIFIC



# AMD 3+

ACCURATE MULTIPLE DELIVERY

## Somni Scientific's AMD 3+ (Accurate Multiple Delivery)

The ultimate in Inhalant Anesthesia gas delivery. This pressurized anesthesia system can deliver accurate adjustable flowrates to 3 stations simultaneously.

Each station is equipped with a flowmeter and toggle switch that allows the user to leave flows preset or adjust to their needs.

Unlike any other system on the market this allows users the greatest flexibility for their procedures.

Any combination of Nose Cones, Induction Chambers, Stereotaxic or Imaging device can be used with the AMD-3+.

### ► SPECIFICATIONS

Weight: 11 lbs / 5 kg

Height: 7 in / 177.8 mm

Depth: 8 in / 203.2 mm

Width: 12 in / 304.8 mm




### FEATURES:

- 3 Flowmeter Controlled Stations with Toggle Switches
- Color Coded Tubing to each station
- Common Gas outlets for each station
- Small footprint
- Adjustable Flowrates
- Stacking design saves valuable workbench space



### INCLUDES:

- TT-2 with Special 6psi Calibrated vaporizer
- 3 each, 6' color coded tubing with adapters



1-877-637-3625

info@somniscientific.com

www.somniscientific.com

Proudly machined, assembled, and calibrated in the USA