

SOMNI  
SCIENTIFIC



## SOMNI REBREATHING HEAD



### SAFETY:

- Pressure Manometer center mount for immediate visibility.
- Pop off valve releases excess pressure at 20cm H2O.
- Large head designed allowing to maintain a safe internal pressure.



### ▶ COMPATIBLE ACCESSORY OPTIONS:

Adult breathing circuit

Pediatric breathing circuit

Pop off Valve

Pop off Valve with occlusion

Pop off occlusion valve




### QUALITY:

- Made with top-of-the-line, durable PCV, with protected powder coating to withstand use.
- High-grade silicone O-rings for a leak proof system.



### CONVENIENT:

- Fitted absorbent canister with our screw top design for ease of use.
- Flutter domes and rebreathing bag in proper view and designed for ease of breath application for user.
- Large fitted absorbent canister fits 3lb bag of Sodalime.



1-877-637-3625

info@somniscientific.com

www.somniscientific.com

Proudly machined, assembled, and calibrated in the USA