

SOMNI

SCIENTIFIC

TT2 ANESTHESIA SYSTEM

TABLE TOP 2* ANESTHESIA SYSTEM WITH BUILT IN FRESH GAS DUAL DELIVERY VALVE




FEATURES:

- Compact 12 inch footprint.
- Compatible with all vaporizer mounts.
- Designed for ease of use, setup, and portability.
- Flowmeter with mechanical stop flow control knob, 0-4 LPM.
- Customizable with the ability to add additional accessories.



SPECIFICATIONS:

- **Weight With Vaporizer:** 25.15 lbs
- **Weight Without Vaporizer:** 9.40 lbs



QUALITY:


- Machined from anodized aluminum for durability and chemical resistance.
- Ability to utilize one or both fresh gas delivery outlets without disconnecting supply lines.
- 5 year warranty.



INCLUDES:

- Vaporizer inlet and outlet endcaps.
- Color coded fresh gas delivery lines with common outlet and 15 mm adaptor ends (2).

**Vaporizer and other accessories not included.*



1-877-637-3625



info@somniscientific.com



www.somniscientific.com

Proudly machined, assembled, and calibrated in the USA