

SOMNI
SCIENTIFIC



THE SOMNI WORK STATION

► SPECIFICATIONS

Weight: 35 lb / 15.88 kg

Height: 40 in / 1016 mm

Depth: 20.5 in / 520.7 mm

Width: 21 in / 533.4 mm



OPTIONAL ACCESSORIES:

- TT2 Anesthesia System
- EPS Active Vacuum System
- AMD-3 Multiple Flowmeter System
- Single or Dual E-tank Block with regulator & hose
- WAG Cannister and Holder
- Iso/Sevo bottle Holder
- Mounted Power Strip
- Magnetic Nose Cones and Manifolds (Miniflow, Uniflow, Multi station manifolds)
- Induction Chamber with Magnetic Mounting

The Somni Work Station

A mobile anesthesia system that can also be utilized as a workspace. It is customizable to suit any researchers needs.

Utilizing an “H” style frame makes the Somni WS very stable and able to fit into tight spaces. With options like E-tanks or an O2 Concentrator mounting system it becomes a fully independent workstation.




*Shown with
Somni W.A.G.,
EPS-3, AMD-
3+, Uniflow
and Induction
Chamber*



STANDARD FEATURES INCLUDE:

- Powder Coated Aluminum Frame, cleans with ease
- Magnetic Work Area for use with Somni Induction Chambers and Nose Cones
- Standing Height work area for Improved Ergonomics
- Shelf for EPS-3 active vacuum system or Cages
- 3” Easy Roll Wheels
- Magnetic Tubing Organizer



1-877-637-3625

info@somniscientific.com

www.somniscientific.com

Proudly machined, assembled, and calibrated in the USA